

Evaluation of an Electron Beam Energy Verification Method Using Statistical Process Control

Stephanie A. Parker, MS, DABR
Novant Health GWSM



Kernersville



Winston-Salem



Salisbury

Conflicts of Interest

- * I have no conflicts of interest to disclose

Outline

- * Background

- * Measurement Methods

- * Specification Limits

- * Statistical Process Control Techniques

 - * Control Limits

 - * Process Capability


 - * Process Acceptability

Background

- * Electron Energy Constancy
 - * Important component of routine linac QA
 - * TG142 – Monthly Check
 - * Challenging Measurement
 - * Rapid falloff of depth dose beyond D-Max
 - * Multiple electron energies per machine

Background

- * Electron Energy Constancy
 - * Measurement method using detector array and wedge shaped filter
 - * Described by several authors
 - * 1991 - 2011
 - * Automated feature of some array detectors



Background

- * Electron Energy Constancy
 - * Equipment may be limited at some facilities
 - * Smaller Clinics
 - * Satellite Facilities
 - * Budget Constraints
 - * Method described applicable to most array detectors

Outline

- * Background
- * Measurement Methods
- * Specification Limits
- * Statistical Process Control Techniques
 - * Control Limits
 - * Process Capability
 - * Process Acceptability

Measurement Method

- * Closely Followed Method Described by Watts in 1998

Evaluation of a diode detector array for use as a linear accelerator QC device

Ronald J. Watts^{a)}

Live Oak Regional Cancer Center, San Antonio, Texas 78233

Med. Phys. 25 (2), February 1998

Measurement Method

- * IBA Blue Phantom
 - * Annual QA
 - * $E_{p,0}$ Measured for Each Electron Beam




Measurement Method

- * MapCheck2 Diode Array
- * Aluminum Wedge

7,07 mm spacing, entire array


1527 total diodes

32,0 x 26,0 cm array size



Power/Data Input

Status Indicators



Inherent Buildup: $2.0 \pm 0.1 \text{ g/cm}^2$

MapCHECK 2

Radiation Measured:

Photons: Co-60 to 25 MV

Electrons: 6 MeV to 25 MeV

Measurement Method

- * Equipment Setup
 - * MapCheck2
 - * Leveled
 - * Centered with CAX
 - * 100 cm SSD
 - * 20 cm x 20 cm Electron Cone




Measurement Method

- * Equipment Setup
 - * Aluminum Wedge
 - * Wedge Direction Inplane
 - * Toe Toward Gantry
 - * Heel Right Angle on MapCheck2 Surface




Measurement Method

- * Equipment Setup
 - * Aluminum Wedge
 - * Toe Aligned with 15 cm Field Edge Demarcation




Measurement Method




Measurement Method

- * Planar Fluences were Measured for Each Electron Energy
- * Only readings along center of detector in Y direction of interest



Measurement Method

- * Intercept of the tangent line to the 50% point on the “toe” end of profile calculated and recorded




Measurement Method

- * 4 Linear Accelerators
 - * 1 Varian 21iX
 - * 2 Varian 21 EX
 - * 1 Varian 21 EX-S
- * 5 Electron Energies per Linac
 - * 6, 9, 12, 16, 20 MeV
- * Total of 20 Electron Beams



Measurement Method



Most Probably Electron Energy at Surface ($E_{p,o}$) vs. Array Detector Intercept

Measurement Method

- * Efficiency
 - * Same setup for each electron beam
 - * Adds about 10 minutes to acquire fluences
 - * Use same setup without wedge to measure profiles
 - * Flatness/Symmetry

Outline

- * Background
- * Measurement Methods
- * Specification Limits
- * Statistical Process Control Techniques
 - * Control Limits
 - * Process Capability
 - * Process Acceptability

Specification Limits

Task Group 142 report: Quality assurance of medical accelerators^{a)}

TABLE II. Monthly.

Procedure	Machine-type tolerance	
	Non-IMRT	IMRT
Dosimetry		
X-ray output constancy		
Electron output constancy		2%
Backup monitor chamber constancy		
Typical dose rate ^a output constancy	NA	2% (@ IMRT dose rate)
Photon beam profile constancy		1%
Electron beam profile constancy		1%
Electron beam energy constancy		2%/2 mm

Specification Limits

- * Need to correlate shift in PDD to change in $E_{p,0}$
- * Referred to TG-70

Recommendations for clinical electron beam dosimetry: Supplement to the recommendations of Task Group 25

Gerbi *et al.*: TG70: Recommendations for clinical electron beam dosimetry

Med. Phys. 36 (7), July 2009

Specification Limits

$E_{p,0}$: Most probable energy (kinetic) of an electron beam at the surface of a water phantom for an electron beam. Unit: MeV.


$$E_{p,0} = 0.22 + 1.98R_p + 0.0025R_p^2$$

$$R_p = 1.271R_{50} - 0.23 \quad (\text{cm}).$$

Need to Know: Change in $E_{p,0}$ for a 2 mm change in R_{50}

Specification Limits

$E_{p,0}$ vs R_{50}



Specification Limits

$$y = 2.5595x - 0.3241$$

- * Slope = 2.56 MeV/cm
- * Slope = 0.256 MeV/mm
- * 2 mm Shift in R_{50} = 0.51 change in MeV
- * Spec = +/- 0.5 MeV

Clinical Implementation

Electron Energy Checks

Setup: 100 cm SSD to Mapcheck Surface, No Buildup, Aluminum Wedge Placed on Mapcheck Surface with Toe toward gantry touching 15 cm line, Difference Between Baseline and Calculated Epo Should be ≤ 0.5 MeV

	6 MeV Electrons		9 MeV Electrons		12 MeV Electrons	
	A	B	A	B	A	B
Detector Location	70	80	60	70	40	50
Readings						
Calculated Intercept	#DIV/0!		#DIV/0!		#DIV/0!	
Calculated Epo	#DIV/0!		#DIV/0!		#DIV/0!	
Baseline	5.769		8.804		11.842	
Difference	#DIV/0!		#DIV/0!		#DIV/0!	

Statistical Process Control Techniques

- * Remaining Questions
 - * Stability of Method Over Time
 - * Reproducibility of Setup
 - * Inter-User Variability
- * Turned to Statistical Process Control Techniques
 - * A lot of options

Statistical Process Control Techniques

- * Most Applicable
 - * Control Limits & Control Charts
 - * Process Capability
 - * Process Acceptability

Control Charts

- * Originated with Walter Shewhart
 - * 1920's
 - * Bell Laboratories
- * Used to determine if Process:
 - * Stable
 - * Has Predictable Performance



Control Charts

- * Shewhart Identified two sources of process variation
 - * “Chance” Variation
 - * Inherent in process
 - * Stable over Time
 - * “Assignable” Variation
 - * Result of specific event outside system
 - * Unstable over Time

Control Charts

- * “Chance” Variation
 - * Random Error
 - * Common Cause
- * “Assignable” Variation
 - * Systematic Error
 - * Special Cause

Control Charts Help Distinguish Between
the Two Types of Error

Control Charts

- * Different Types of Control Charts
 - * Attribute Data
 - * Discrete
 - * Y/N
 - * Good/Bad
 - * Variable Data
 - * Continuous Scale

Control Charts

- * Variable Data Charts (Actually Pairs of Charts)

- * X and Moving Range Chart

- * Sample Size (n) = 1


- * X-Bar and Range Chart

- * n = 2-9

- * X-Bar and S Chart

- * n > 10

Control Charts



Control Charts

- * How do I apply control charts to $E_{p,0}$ measurements?

Statistical process control for radiotherapy quality assurance

Todd Pawlicki^{a)}

Department of Radiation Oncology, Stanford University School of Medicine, Stanford, California 94305

Matthew Whitaker

Radiological Imaging Technology, Inc., 637 Elkton Drive, Colorado Springs, Colorado 80907

Arthur L. Boyer

Department of Radiation Oncology, Stanford University School of Medicine, Stanford, California 94305

Med. Phys. 32 (9), September 2005

Control Charts

- * Steps
 - * Collect Initial Data
 - * Ensure meets specification
 - * Establish Control Limits Using Collected Data
 - * Create Charts
 - * Record Data in Charts

Control Charts

- * 6 months of data collected for each electron beam
 - * Calculated the Average (\bar{x})
 - * Calculated Range (R)
 - * Absolute difference between two successive measurements
 - * Calculated Average Range (\bar{R})

Control Charts

- * *Control Limits for Individual Control Chart*

- * *Control Center*(C_c) = \bar{x}

- * *Upper Control Limit*(UCL) = $\bar{x} + 2.66\bar{R}$

- * *Lower Control Limit* (LCL) = $\bar{x} - 2.66\bar{R}$

- * *Control Limits for Moving Range Chart*

- * *Range Center*(R_c) = \bar{R}


- * *Upper Range Limit*(URL) = $3.27\bar{R}$

- * *Lower Range Limit*(LRL) = 0




**Control Limits
Represent 3
Standard Errors
from the Mean**

Individual Control Chart




Moving Range Chart




Control Charts

- * All control limits well within specification limits
- * Shows that the process is in control
 - * All data points fall within control limits
 - * Data follows a random pattern
- * Process stability should allow for distinction between random and systematic errors
 - * Noise small so signal should be able to be detected

Individual Control Chart



Moving Range Chart



Statistical Process Control Techniques

- * Most Applicable
 - * Control Limits & Control Charts
 - * Process Capability
 - * Process Acceptability

Process Capability and Acceptability

JOURNAL OF APPLIED CLINICAL MEDICAL PHYSICS, VOLUME 14, NUMBER 1, 2013

Retrospective analysis of linear accelerator output constancy checks using process control techniques


Taweap Sanghangthum,^{1,2} Sivalee Suriyapee,² Somyot Srisatit,¹
Todd Pawlicki^{3,a}

*Department of Nuclear Technology,¹ Chulalongkorn University, Bangkok, Thailand;
Division of Radiation Oncology,² King Chulalongkorn Memorial Hospital, Bangkok,
Thailand; Department of Radiation Medicine and Applied Sciences,³ University of
California, San Diego, USA
tpaw@ucsd.edu*

Process Capability and Acceptability

- * Process Capability (C_p)
 - * How well a process is capable of meeting specification
 - * Comparison between:
 - * Spread of data
 - * Window size of specification limits
 - * $C_p = \frac{USL - LSL}{6\sigma}$
 - * σ = process standard deviation

Process Capability and Acceptability




- * Process Capability (C_p)
- * $C_p = \frac{USL - LSL}{6\sigma}$
- * If $C_p \geq 1$, process capable of meeting specs

Process Capability and Acceptability


- * Process Acceptability (C_{pk})
 - * How well a process is centered within the specification limits
 - * $C_{pk} = \min\left(\frac{USL - \mu}{3\sigma}, \frac{\mu - LSL}{3\sigma}\right)$
 - * μ = process mean
 - * σ = process standard deviation

Process Capability and Acceptability



- * Process Acceptability (C_{pk})
 - * $C_{pk} = \min\left(\frac{USL - \mu}{3\sigma}, \frac{\mu - LSL}{3\sigma}\right)$
 - * If $C_{pk} \geq 1$, process centered within specs
 - * If process perfectly centered, two solutions would be equal

Process Capability and Acceptability




Sanghangthum et al, JACMP, Vol 14 (1),
2013, Pg 151, Used with Permission

Process Capability and Acceptability

- * Data Requirements
 - * Normal Distribution
 - * “Large Enough” data sample
 - * Typically ≥ 30 data points
 - * Need to assess measured data to determine if requirements met

Process Capability and Acceptability

- * Normal Distribution Assessment
 - * Anderson-Darling Test
 - * $P \geq 0.05$, data considered normal
 - * P values range 0.14 – 1.0
 - * All data considered normally distributed
- * Data Size ~ 40



Conclusion: Data Meets Requirements for Analysis Using C_p and C_{pk}

Process Capability and Acceptability

- * C_p Results
 - * All values > 1
 - * 2.7 – 10
 - * Process for all energies capable of meeting specifications

Process Capability and Acceptability

- * C_{pk} Results

- * All values > 1


- * 1.8 – 9.6

- * Center of the process for all energies within specification limits

Process Capability and Acceptability

- * Comparison of C_p and C_{pk}
 - * $C_{pk} < C_p$ in all but 1 case
 - * Indicated that process has some shift from baseline
 - * Baseline value single measurement
 - * Indicates that target value should be an average of in control values vs. single measurement

Process Capability and Acceptability



Summary

- * Measurement Method Established to Measure $E_{p,0}$
 - * MapCheck2 and Aluminum Wedge
- * Specification Limits Established
 - * 2 mm PDD shift = +/- 0.5 MeV
- * Statistical Process Control Techniques
 - * Control Limits
 - * Process Capability
 - * Process Acceptability

Summary

- * Statistical Process Control Techniques Utilized
 - * Control Limits
 - * Process in Control
 - * Control Limits well within Specification Limits
 - * Process Capability
 - * $C_p > 1$ in all Cases
 - * Process Capable of Meeting Specs
 - * Process Acceptability
 - * $C_{pk} > 1$ in all Cases
 - * Process Centered within Specification Limits


Conclusion

- * Electron energy verification method
 - * Efficient
 - * Effective
- * Good option for centers with limited equipment
 - * Small Centers
 - * Satellite Centers
- * Statistical Process Control
 - * Tools useful for analyzing QA processes

La Fin

Merci

Me at the Marie Curie
Museum in Paris



Questions?



References

- * R. J. Watts, "Evaluation of a diode detector array for use as a linear accelerator QC device," *Med Phys* 25, 247-250 (1998).
- * R. P. King and R.S. Anderson, "A simple method for electron energy constancy measurement," *J. Appl. Clin. Med. Phys.* 2, 51-53 (2001).
- * U. F. Rosenow, M. K. Islam, H. Gaballa, and H. Rashid, "Energy constancy checking for electron beams using a wedge-shaped solid phantom combined with a beam profile scanner," *Med. Phys.* 18, 19-25, (1991).
- * B. J. Gerbi, J. A. Antolak, F. C. Deibel, D. S. Followill, M. G. Herman, P. D. Higgins, M. S. Huq, D. N. Mihalidis, and E. D. Yorke, "Recommendations for clinical electron beam dosimetry: Supplement to the recommendations of Task Group 25," *Med. Phys.* 36, 3239-3279 (2009).
- * E. E. Klein, J. Hanley, J. Bayouth, F. Yin, W. Simon, S. Dresser, C. Serago, F. Aguirre, L. Ma, B. Arjomandy, and C. Liu, "Task Group 142 report: Quality assurance of medical accelerators," *Med. Phys.* 36, 4197-4212 (2009).
- * T. Pawlicki, M. Whitaker, A. L. Boyer, "Statistical process control for radiotherapy quality assurance," *Med. Phys.* 32, 2777-2786 (2005).
- * T. Sanghangthum, S. Suriyapee, S. Srisatit, T. Pawlicki, "Retrospective analysis of linear accelerator output constancy checks using process control techniques," *J. Appl. Clin. Med. Phys.* 14, 147-160 (2013).
- * Montgomery, Douglas, *Introduction to Statistical Process Control*, 5th Edition, 2005, ISBN 0-471-65631-3
- * NIST/SEMATECH e-Handbook of Statistical Methods, <http://www.itl.nist.gov/div898/handbook/>